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Viscosity and Viscoelastic Characterization

From early stage formulation to product use

Ease of use next generation rheometer
Innovative microfluidic rheometry
process-relevant material testing
Industry standard Asphalt binder grade
low volume viscosity and size analysis
Relative dilute polymer solutions viscosity

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Bio/Pharmaceuticals

Antibody Therapeutics
 Drug Injectability
 Protein Stability
 Protein Therapeutics
 RNA Therapeutics
 Eye drops
 Blood Plasma

Rechargeable Batteries

Lithium Perchlorate electrolyte
 PVDF/HP solution in NMP

Conductive/Ceramic ink

Inkjet ink_ high Shear Rate Viscosity

Food

Infant Milk

Oil Lubricant

Cosmetics/Detergent/Nano Particles

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
Method (# of Samples)	microVISC™	Glass Capillary
Multi-point (6)	25	120
Single point (2)	6	20

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Measurement Protocol:

Loading...

Load your sample into one of *microVISC*® disposable pipettes. Place it in the system and you are ready!

Simple!

Measuring...

Perform a quick automatic measurement or run at several shear rates using advanced mode.

Precise!

Cleaning...

After testing or between insoluble samples run a pipette of appropriate solvent using cleaning mode.

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

- Verify the quality of your oil quickly
Minimum sample: 100 µL
- Humidity Range 0 - 90%, non-condensing
- Viscosity Range 1- 3500 cSt
- Accuracy +/- 2 % Reading
Repeatability +/- 0.5% Reading
- Operating Range 18 - 40 °C
Output Range 20 - 100 °C
- <1 minute
Fast viscosity measurements

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K = ball constant mPa·s·cm³/g·s

ρ_1 = density of the ball in g/cm³

ρ_2 = density of the liquid to be measured at the measuring temperature in g/cm³


t = falling time of the ball in seconds.

Ball N ^o	Viscosity range (mPa·s)	Material	Density (g/cm ³)
1	0.6 to 10	Borosilicate glass	2.2
2	7 to 130	Borosilicate glass	2.2
3	30 to 700	Nickel iron alloy	8.1
4	200 to 4,800	Nickel iron alloy	8.1
5	1,500 to 45,000	Stainless steel	7.7-8.1
6	> 7,500	Stainless steel	7.7-8.1
G	Gas Viscosity	Borosilicate glass	2.2

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ISO 2431

Code	Bore Ø mm	Viscosity range (cSt)
FCI003	3	5 to 42
FCI004	4	35 to 135
FCI005	5	100 to 350
FCI006	6	190 to 680

ASTM D-1200

Code	Bore Ø mm	Viscosity range (cSt)
FCA002	2.53	25 to 120
FCA003	3.40	40 to 220
FCA004	4.12	70 to 370
FCA005	5.20	200 to 1200

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
AASHTO 規範：

- Pass/Fail medium temp. original binder test T315-11)
- Pass/Fail high temperature RTFO binder (T315-11)
- Pass/Fail high temperature PAV binder (T315-11)
- Linearity test (T315-12)

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- 鋪面等級判斷 (Research Grade Software)
- Grade determination test (R29-02)

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

- 3M
- BASF
- Eastman
- US Army Corp of Engineers
- Texas A&M University
- Georgia Institute of Technology
- Green Mountain Coffee NIST
- Bausch & Lomb

最小分析黏度 cP	1cP (mPa.s)*
最高分析黏度 cP	100M cP (mPa.s)*
轉速範圍 RPM	0.1 to 2000 1/min
角速度範圍 Angular Velocity	0.01 to 200 rad/s*
剪切速率範圍 Shear Rate	Measuring System Dependent
剪切力範圍 Shear Stress	Measuring System Dependent
扭力範圍 Torque	0.001 to 20 mNm
內建溫控精度 °C	0.01°C
內建溫控範圍 °C	-10 to 120 °C
溫度控制間隔	+/- 0.1 °C
夾具 Geometry	Cone & Plate, Parallel, Plate, and Co-Axial Cylinders...

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Thermostatic Bath

PPR Plugging System

Universal adapter for all Brookfield

viscometer

	ALPHA	SMART	EXPERT	PREMIUM
Fluid	Recommended for Newtonian fluids	All types	All types	All types
Description	Fast and accurate viscosity readings. Easy to use and low budget	Allows temperature control and data register. Indispensable on QC and R&D laboratories	Add essential performances to determine viscosity and other rheological features	Viscometer monitored by our DataBoss software, offers a wide range and unique range of rheological applications
LCD Display	✓ alphanumeric	✓ graphic display	✓ graphic display	✓ graphic display
Touch keyboard	✗	6 keys touch keyboard	12 keys touch keyboard	12 keys touch keyboard
Temperature Probe PT100	✗	Optional	✓	✓
Temperature Settings	✗	°C/°F (optional)	°C/°F	°C/°F
Datalogger Software	✗	✓	✓	✓
DataBoss Software	✗	✗	✗	Under request
Bidirectional Unit	✗	✗	✗	✓
10 Languages	✓	✓	✓	✓
Working memories	✗	9	9	9
Speeds	18	18	54	2600
Speed settings	✗	✗	Customizable by user	Customizable by user
RPM	0.3-100	0.3-100	0.01-200	0.01-250
Data displayed	Speed Spindle Viscosity (cP or mPa.s) % of torque	Speed Spindle Viscosity (cP or mPa.s) % of torque Shear rate (w.coaxial spindles)(s ⁻¹) Shear stress (with coaxial spindles)(N/m ²) Density (g/cm ³)	Speed Spindle Viscosity (cP or mPa.s) % of torque Shear rate (w.coaxial spindles)(s ⁻¹) Shear stress (with coaxial spindles)(N/m ²) Density (g/cm ³)	S Speed Spindle Viscosity (cP or mPa.s) % of torque Shear rate (w.coaxial spindles)(s ⁻¹) Shear stress (with coaxial spindles)(N/m ²) Density (g/cm ³)
Density cSt	✗	By default and introduced by user	By default and introduced by user	By default and introduced by user
Viscosity readings	Absolute viscosity: cP or mPa.s Apparent viscosity: cP or mPa.s	Absolute viscosity: cP or mPa.s Apparent viscosity: cP or mPa.s Kinematic viscosity cSt r mm ² /sec	Absolute viscosity: cP or mPa.s Apparent viscosity: cP or mPa.s Kinematic viscosity cSt r mm ² /sec	Absolute viscosity: cP or mPa.s Apparent viscosity: cP or mPa.s Kinematic viscosity cSt r mm ² /sec
Time Settings	✗	✓	✓	✓
Program settings	✗	Time to torque/Time to stop	Time to torque/Time to stop Step/ Multistep/ Ramp	Time to torque/Time to stop Step/ Multistep/ Ramp
Competing products	DVE and Dial reading	DVI Prime	DVII Pro	DVIII Ultra

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